

Rabbit Anti-Human IL9 Monoclonal Antibody For If/Icc And Fc (Intra) - P15248

Item Number: CM0006234



Introduction

-conjugated rabbit anti-human IL9 monoclonal antibody validated for immunofluorescence (IF/ICC) and intracellular flow cytometry (FC). Targets interleukin-9, a cytokine involved in immune response and inflammation. For human samples. SWISS: P15248.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL9 mAb
Remarks/Alias	IL9; HP40; IL-9; P40; interleukin-9
Species	Human
GeneID (Human)	3578
GeneID	3578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-144 of human IL9(NP_000581.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC (intra)
Cross Reactivity	Human
SWISS	P15248
Protein Weight	15kDa
Shipping	Ice bag