

Rabbit Anti-Human IL9 Monoclonal Antibody For If/Icc, Fc (Intra) - P15248

Item Number: CM0005759



Introduction

-conjugated rabbit monoclonal antibody against human IL9, validated for IF/ICC and intracellular flow cytometry. Ideal for studying IL9-mediated immune responses, T-cell differentiation, and mast cell biology.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL9 mAb
Remarks/Alias	IL9; HP40; IL-9; P40; interleukin-9
Species	Human
Gene ID (Human)	3578
Immunogen	Recombinant fusion protein containing amino acids 19-144 of human IL9 (NP_000581.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC (intra)
Cross Reactivity	Human
SWISS	P15248
Protein Weight	15kDa
Shipping	Ice bag