

Rabbit Anti-Human Integrin Alpha 3/Cd49C Mab - P26006

Item Number: CM0010511

Introduction

rabbit anti-human Integrin alpha-3 (CD49c) mAb for FC. Targets the laminin/collagen receptor; validated for human samples. Ideal for cell adhesion, migration, and invasion studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Integrin alpha 3/CD49c mAb
Remarks/Alias	JEB7; VL3A; CD49C; FRP-2; GAPB3; ILNEB; MSK18; VCA-2; VLA3a; GAP-B3
Species	Human
GeneID (Human)	3675
GeneID	3675
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 466-620 of human Integrin alpha 3/CD49c (NP_002195.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26006
Protein Weight	117kDa
Shipping	Ice bag