

Rabbit Anti-Human Integrin-B1/Cd29 Monoclonal Antibody For Fc - P05556

Item Number: CM0004296



Introduction

-conjugated rabbit monoclonal antibody targeting human Integrin beta-1 (CD29) for flow cytometry. Recognizes human CD29, molecular weight ~87-91 kDa. Ideal for cell adhesion and integrin signaling research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Integrin-β1/CD29 mAb
Remarks/Alias	CD29; FNRB; MDF2; VLAB; GPIIA; MSK12; VLA-BETA
Species	Human
Gene ID	3688
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-728 of human Integrin-β1/CD29 (NP_002202.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P05556
Protein Weight	87kDa/88kDa/91kDa
Shipping	Ice bag