


Rabbit Anti-Human Integrin B5/Itgb5 Monoclonal Antibody For Flow Cytometry - P18084

Item Number: CM0009783



Introduction

-conjugated rabbit monoclonal antibody targeting human Integrin $\beta 5$ (ITGB5). Validated for flow cytometry, this antibody recognizes the fibronectin receptor, suitable for cell adhesion studies and viral entry research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Integrin $\beta 5$ /ITGB5 mAb
Remarks/Alias	ITGB5; integrin beta-5
Species	Human
GeneID	3693
Immunogen	Recombinant protein corresponding to a sequence within amino acids 24-490 of human Integrin $\beta 5$ /ITGB5 (NP_002204.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P18084
Protein Weight	88kDa
Shipping	Ice bag