


Rabbit Anti-Human Itga1/Cd49A Mab For Fc - P56199

Item Number: CM0010597



Introduction

High-quality Rabbit anti-Human ITGA1/CD49a monoclonal antibody for flow cytometry. Validated for human samples, recognizes integrin alpha-1 (CD49a/VLA-1), ideal for cell adhesion and immune cell studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ITGA1/CD49a mAb
Remarks/Alias	VLA1; CD49a
Species	Human
GeneID (Human)	3672
GeneID	3672
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-1141 of human ITGA1/CD49a (NP_852478.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P56199
Protein Weight	131kDa
Shipping	Ice bag