


## Rabbit Anti-Human Itga1/Cd49A Monoclonal Antibody For Fc - P56199

Item Number: CM0008379



### Introduction

-conjugated rabbit monoclonal antibody against human integrin alpha-1 (ITGA1/CD49a/VLA-1), validated for flow cytometry (FC). Suitable for studying cell adhesion and integrin signaling pathways.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human ITGA1/CD49a mAb                                                                                        |
| Remarks/Alias    | VLA1; CD49a                                                                                                              |
| Species          | Human                                                                                                                    |
| GeneID (Human)   | 3672                                                                                                                     |
| GeneID           | 3672                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 29-1141 of human ITGA1/CD49a (NP_852478.1) |
| Source           | Rabbit                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                    |
| Application      | FC                                                                                                                       |
| Cross Reactivity | Human                                                                                                                    |
| SWISS            | P56199                                                                                                                   |
| Protein Weight   | 131kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                  |