

Rabbit Anti-Human Jam-A/Cd321/F11R Monoclonal Antibody For Flow Cytometry - Q9Y624

Item Number: CM0010796



Introduction

-conjugated rabbit monoclonal antibody targeting human JAM-A (CD321, F11R) validated for flow cytometry. Suitable for studying tight junction dynamics, immune cell transmigration, and platelet activation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human JAM-A/CD321/F11R mAb
Remarks/Alias	JAM, KAT, JAM1, JAMA, JCAM, CD321, PAM-1
Species	Human
Gene ID	50848
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-242 of human JAM-A/CD321/F11R (NP_058642.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Human
SWISS	Q9Y624
Protein Weight	27kDa/32kDa
Shipping	Ice bag