

Rabbit Anti-Human Jam-A/Cd321/F11R Monoclonal Antibody For Flow Cytometry - Q9Y624

Item Number: CM0009247



Introduction

Rabbit anti-Human JAM-A/CD321/F11R monoclonal antibody for flow cytometry. Recognizes human JAM-A, a tight junction protein involved in leukocyte transmigration. Recombinant protein immunogen, suitable for FC. For research use.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human JAM-A/CD321/F11R mAb
Remarks/Alias	JAM; KAT; JAM1; JAMA; JCAM; CD321; PAM-1
Species	Human
GeneID	50848
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-242 of human JAM-A/CD321/F11R (NP_058642.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9Y624
Protein Weight	27kDa/32kDa
Shipping	Ice bag