


Rabbit Anti-Human Kir2DL2/Kir2DL3 Monoclonal Antibody For If/Icc And Fc - P43628

Item Number: CM0012446



Introduction

-conjugated rabbit monoclonal antibody specific for human KIR2DL2/KIR2DL3, a killer-cell immunoglobulin-like receptor. Validated for IF/ICC and FC. Ideal for studying NK cell inhibitory signaling and HLA-C recognition.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human KIR2DL2/KIR2DL3 mAb
Remarks/Alias	KIR2DL3; CD158B2; CD158b; GL183; KIR-023GB; KIR-K7b; KIR-K7c; KIR2DS5; KIRCL23; NKAT; NKAT2; NKAT2A; NKAT2B; p58; killer cell immunoglobulin-like receptor 2DL3
Species	Human
GeneID (Human)	3803/3804
GeneID	3803/3804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-245 of human KIR2DL3 (NP_056952.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	P43628
Protein Weight	27kDa/37kDa/38kDa
Shipping	Ice bag