


Rabbit Anti-Human Kir2DL2/Kir2DL3 Monoclonal Antibody For If/Icc, Fc - P43628

Item Number: CM0005078



Introduction

High-performance -conjugated rabbit monoclonal antibody against human KIR2DL2/KIR2DL3 (P43628). Suitable for immunofluorescence and flow cytometry. Validated for NK cell immunophenotyping and inhibitory receptor studies.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human KIR2DL2/KIR2DL3 mAb
Remarks/Alias	KIR2DL3; CD158B2; CD158b; GL183; KIR-023GB; KIR-K7b; KIR-K7c; KIR2DS5; KIRCL23; NKAT; NKAT2; NKAT2A; NKAT2B; p58; killer cell immunoglobulin-like receptor 2DL3
Species	Human
GeneID (Human)	3803/3804
GeneID	3803/3804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-245 of human KIR2DL3(NP_056952.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	P43628
Protein Weight	27kDa/37kDa/38kDa
Shipping	Ice bag