

Rabbit Anti-Human Kir2DL4/Cd158D Mab For Fc - Q99706

Item Number: CM0012618



Introduction

-conjugated rabbit monoclonal antibody targeting human KIR2DL4/CD158d. Validated for flow cytometry (FC). Recognizes extracellular domain (aa 24-242). Suitable for NK cell receptor studies and immunology research.

[Learn More](#)

Parameter	Details
Product Name	Rabbit anti-Human KIR2DL4/CD158d mAb
Remarks/Alias	KIR2DL4; CD158D; G9P; KIR-103AS; KIR-2DL4; KIR103; KIR103AS; killer cell immunoglobulin-like receptor 2DL4
Species	Human
GeneID (Human)	3805
GeneID	3805
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-242 of human KIR2DL4/CD158d (NP_002246.5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q99706
Protein Weight	24kDa/33kDa/35kDa/37kDa/39kDa/41kDa
Shipping	Ice bag