

Rabbit Anti-Human Kir2DL4/Cd158D Mab For Fc - Q99706

Item Number: CM0004987



Introduction

Rabbit monoclonal anti-human KIR2DL4/CD158d antibody conjugated to . Validated for flow cytometry, targeting the killer cell immunoglobulin-like receptor that recognizes HLA-G and mediates NK cell functions at the maternal-fetal interface.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human KIR2DL4/CD158d mAb
Remarks/Alias	KIR2DL4; CD158D; G9P; KIR-103AS; KIR-2DL4; KIR103; KIR103AS; killer cell immunoglobulin-like receptor 2DL4
Species	Human
GeneID (human)	3805
GeneID	3805
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-242 of human KIR2DL4/CD158d (NP_002246.5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q99706
Protein Weight	24kDa/33kDa/35kDa/37kDa/39kDa/41kDa
Shipping	Ice bag