

Rabbit Anti-Human Klrg1 Polymab® For Fc - Q96E93

Item Number: CM0005094



Introduction

-conjugated rabbit anti-human KLRG1 monoclonal antibody for flow cytometry. Reacts specifically with human KLRG1, an inhibitory receptor on NK and T cells. Ideal for immune cell profiling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human KLRG1 PolymAb®
Remarks/Alias	2F1; MAFA; MAFA-L; CLEC15A; MAFA-2F1; MAFA-LIKE
Species	Human
GeneID	10219
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-195 of human KLRG1 (NP_001316028.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q96E93
Protein Weight	22kDa
Shipping	Ice bag