

Rabbit Anti-Human Klrp1 Polymab® For Flow Cytometry - Q96E93

Item Number: CM0009812



Introduction

-conjugated recombinant rabbit monoclonal antibody targeting human KLRG1 for flow cytometry. Recognizes a 22 kDa inhibitory receptor expressed on NK cells and CD8+ T-cells. Validated in human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human KLRG1 PolymAb®
Remarks/Alias	2F1; MAFA; MAFA-L; CLEC15A; MAFA-2F1; MAFA-LIKE
Species	Human
Gene ID	10219
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-195 of human KLRG1 (NP_001316028.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q96E93
Protein Weight	22kDa
Shipping	Ice bag