


Rabbit Anti-Human Lambda Light Chain Monoclonal Antibody For Flow Cytometry - P0CG04

Item Number: CM0008083



Introduction

-conjugated rabbit monoclonal antibody targeting human immunoglobulin lambda light chain constant region. Validated for flow cytometry (FC) analysis of human samples. Ideal for B cell and humoral immunity research; recombinant immunogen (aa 1-106), ~11kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Lambda Light Chain mAb
Species	Human
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-106 of human Lambda Light Chain
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P0CG04
Protein Weight	11kDa
Shipping	Ice bag