


Rabbit Anti-Human Lambda Light Chain Monoclonal Antibody For Flow Cytometry - P0Cg04

Item Number: CM0005319



Introduction

High-quality rabbit monoclonal antibody against human lambda light chain, conjugated to for flow cytometry. Validated for human samples in immunological research. Swiss-Prot: P0CG04.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human Lambda Light Chain mAb
Species	Human
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-106 of human Lambda Light Chain
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P0CG04
Protein Weight	11kDa
Shipping	Ice bag