


Rabbit Anti-Human Lamp1/Cd107A Mab - P11279

Item Number: CM0009494



Introduction

-conjugated rabbit monoclonal antibody against human LAMP1/CD107a for intracellular flow cytometry. Validated for human samples, targets lysosome-associated membrane glycoprotein 1. Suitable for lysosome biogenesis, autophagy, and NK cell cytotoxicity research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human LAMP1/CD107a mAb
Remarks/Alias	LAMPA; CD107a; LGP120
Species	Human
GeneID (Human)	3916
GeneID	3916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-382 of human LAMP1/CD107a (NP_005552.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P11279
Protein Weight	38kDa/45kDa
Shipping	Ice bag