


Rabbit Anti-Human Lamp1/Cd107A Monoclonal Antibody For Flow Cytometry (Intracellular) - P11279

Item Number: CM0008904



Introduction

-conjugated rabbit monoclonal antibody targeting human LAMP1/CD107a. Designed for intracellular flow cytometry (FC intra) to detect a 38–45 kDa lysosomal membrane glycoprotein involved in lysosome biogenesis and NK cell cytotoxicity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human LAMP1/CD107a mAb
Remarks/Alias	LAMPA; CD107a; LGP120
Species	Human
Gene ID (Human)	3916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-382 of human LAMP1/CD107a (NP_005552.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P11279
Protein Weight	38kDa/45kDa
Shipping	Ice bag