


Rabbit Anti-Human Lamp1/Cd107A Monoclonal Antibody For Intracellular Flow Cytometry - P11279

Item Number: CM0008244



Introduction

-conjugated rabbit monoclonal antibody targeting human LAMP1/CD107a (SWISS: P11279) for intracellular flow cytometry. Ideal for lysosomal function, autophagy, and NK cell cytotoxicity research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human LAMP1/CD107a mAb
Remarks/Alias	LAMPA; CD107a; LGP120
Species	Human
GeneID	3916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-382 of human LAMP1/CD107a (NP_005552.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P11279
Protein Weight	38kDa/45kDa
Shipping	Ice bag