

Rabbit Anti-Human Lamp2 Mab For If/Icc And Fc (Intra) - P13473

Item Number: CM0004520



Introduction

-conjugated rabbit anti-human LAMP2 monoclonal antibody, validated for IF/ICC and intracellular flow cytometry. Targets lysosomal membrane protein LAMP2 (CD107b) for autophagy, lysosome biogenesis, and lysosomal storage disease research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human LAMP2 mAb
Remarks/Alias	DND, LAMPB, CD107b, LAMP-2, LGP-96, LGP110
Species	Human
GeneID (Human)	3920
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-375 of human LAMP2 (NP_002285.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC (intra)
Cross Reactivity	Human
SWISS	P13473
Protein Weight	44kDa/45kDa
Shipping	Ice bag