

Rabbit Anti-Human Lamp2 Monoclonal Antibody For If/Icc And Intracellular Fc - P13473

Item Number: CM0011700



Introduction

-conjugated rabbit monoclonal antibody against human LAMP2 (CD107b). Validated for IF/ICC and intracellular flow cytometry. Reacts with human samples. Targets lysosome-associated membrane glycoprotein 2, involved in autophagy and lysosomal function.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human LAMP2 mAb
Remarks/Alias	DND; LAMPB; CD107b; LAMP-2; LGP-96; LGP110
Species	Human
GeneID (Human)	3920
GeneID	3920
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC (intra)
Cross Reactivity	Human
SWISS	P13473
Protein Weight	44kDa/45kDa
Shipping	Ice bag