

Rabbit Anti-Human Lamp2 Monoclonal Antibody For Intracellular Flow Cytometry - P13473

Item Number: CM0010713



Introduction

-conjugated rabbit monoclonal antibody targeting human LAMP2, validated for intracellular flow cytometry. Ideal for lysosome biogenesis, autophagy, and chaperone-mediated autophagy studies. Species reactivity: human. SWISS: P13473.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human LAMP2 mAb
Remarks/Alias	DND; LAMPB; CD107b; LAMP-2; LGP-96; LGP110
Species	Human
GeneID (Human)	3920
GeneID	3920
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-375 of human LAMP2 (NP_002285.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P13473
Protein Weight	44kDa/45kDa
Shipping	Ice bag