

Rabbit Anti-Human Lifr/Cd118 Monoclonal Antibody For Fc - P42702

Item Number: CM0010701



Introduction

Rabbit monoclonal antibody against human LIFR/CD118, directly conjugated with ABflo 647 for flow cytometry (FC). Recognizes human LIFR, a ~124 kDa cell membrane receptor involved in cytokine signaling. Suitable for research use.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human LIFR/CD118 mAb
Remarks/Alias	SWS; SJS2; STWS; CD118; LIF-R
Species	Human
Gene ID (Human)	3977
Gene ID	3977
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P42702
Protein Weight	124kDa
Shipping	Ice bag