

Rabbit Anti-Human LILRA2/Cd85H Monoclonal Antibody For IF/ICC, Fc - Q8N149

Item Number: CM0008881



Introduction

-conjugated rabbit monoclonal antibody targeting human LILRA2/CD85h (ILT1/LIR7). Validated for IF/ICC and flow cytometry. Reacts with human. Ideal for innate immunity and microbial infection research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human LILRA2/CD85h mAb
Remarks/Alias	ILT1; LIR7; CD85H; LIR-7
Species	Human
Gene ID	11027 (Human)
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	Q8N149
Protein Weight	53kDa/51kDa/48kDa/47kDa
Shipping	Ice bag