

Rabbit Anti-Human Lilrb4/Cd85K/Ilt3 Monoclonal Antibody For Fc - Q8Nhj6

Item Number: CM0012707



Introduction

-conjugated rabbit monoclonal antibody against human LILRB4/CD85k/ILT3 (Q8NHJ6). Designed for flow cytometry, targets an inhibitory immune receptor involved in immune tolerance and monocyte regulation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human LILRB4/CD85k/ILT3 mAb
Remarks/Alias	B4, ILT3, LIR5, CD85K, ILT-3, LIR-5
Species	Human
GeneID (human)	11006
GeneID	11006
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8NHJ6
Protein Weight	49kDa
Shipping	Ice bag