


Rabbit Anti-Human Megalin/Lrp2 Monoclonal Antibody For Flow Cytometry (Intracellular) - P98164

Item Number: CM0008838



Introduction

High-quality rabbit monoclonal antibody against human Megalin/LRP2, conjugated with . Validated for intracellular flow cytometry to study this large endocytic receptor in kidney, placenta, and brain. SWISS: P98164.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Megalin/LRP2 mAb
Remarks/Alias	LRP2; DBS; GP330; LDL receptor related protein 2
Species	Human
GeneID (Human)	4036
GeneID	4036
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 4446-4655 of human Megalin/LRP2 (NP_004516.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P98164
Protein Weight	521kDa
Shipping	Ice bag