


Rabbit Anti-Human Megalin/Lrp2 Monoclonal Antibody For If-P And Fc (Intra) - P98164

Item Number: CM0007774



Introduction

High-quality rabbit monoclonal antibody against human Megalin/LRP2, conjugated to . Validated for immunofluorescence and intracellular flow cytometry. Cross-reacts with human and mouse. Ideal for endocytic receptor studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Megalin/LRP2 mAb
Remarks/Alias	LRP2; DBS; GP330; LDL receptor related protein 2
Species	Human
GeneID (Human)	4036
GeneID	4036
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 4446-4655 of human Megalin/LRP2 (NP_004516.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, FC (intra)
Cross Reactivity	Human,Mouse
SWISS	P98164
Protein Weight	521kDa
Shipping	Ice bag