


## Rabbit Anti-Human MERTK Mab For Fc - Q12866

Item Number: CM0008095



### Introduction

-conjugated Rabbit anti-Human MERTK monoclonal antibody validated for flow cytometry (FC). Reacts with Human MERTK. Recombinant protein immunogen. Ideal for apoptosis, phagocytosis, and kinase signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human MERTK mAb                                                                                       |
| Remarks/Alias    | MERTK; MER; RP38; Tyro12; c-Eyk; c-mer; tyrosine-protein kinase Mer                                               |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 10461                                                                                                             |
| Gene ID          | 10461                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 26-499 of human MERTK (NP_006334.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | FC                                                                                                                |
| Cross Reactivity | Human                                                                                                             |
| SWISS            | Q12866                                                                                                            |
| Protein Weight   | 110kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                           |