


Rabbit Anti-Human MERTK Monoclonal Antibody For Fc - Q12866

Item Number: CM0008044



Introduction

-conjugated rabbit monoclonal antibody against human MERTK (MerTK), a receptor tyrosine kinase. Validated for flow cytometry. Targets cell survival, migration, and phagocytosis pathways. Ideal for human cancer and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human MERTK mAb
Remarks/Alias	MERTK; MER; RP38; Tyro12; c-Eyk; c-mer; tyrosine-protein kinase Mer
Species	Human
GeneID (Human)	10461
GeneID	10461
Immunogen	Recombinant fusion protein containing amino acids 26-499 of human MERTK (NP_006334.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q12866
Protein Weight	110kDa
Shipping	Ice bag