

Rabbit Anti-Human Mica/Micb Mab - Q29983 / Q29980

Item Number: CM0004819



Introduction

-conjugated rabbit monoclonal antibody against human MICA/MICB, suitable for flow cytometry. Targets stress-induced NKG2D ligands expressed on various epithelial and tumor cells, aiding immune surveillance studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human MICA/MICB mAb
Remarks/Alias	MIC-A; PERB11.1; PERB11.2; MICB
Species	Human
Gene ID	100507436/4277
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-298 of human MICB (NP_005922.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q29983/Q29980
Protein Weight	4kDa/31kDa/37kDa/43kDa
Shipping	Ice bag