


Rabbit Anti-Human Mmp9 Monoclonal Antibody For Flow Cytometry (Intracellular) - P14780

Item Number: CM0004099



Introduction

-conjugated rabbit monoclonal antibody against human MMP9 for intracellular flow cytometry. Validated for human samples, detects the 78kDa gelatinase B protein involved in extracellular matrix remodeling and leukocyte migration.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human MMP9 mAb
Remarks/Alias	MMP9; CLG4B; GELB; MANDP2; MMP-9; matrix metallopeptidase 9
Species	Human
GeneID (Human)	4318
GeneID	4318
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-707 of human MMP9 (NP_004985.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P14780
Protein Weight	78kDa
Shipping	Ice bag