

Rabbit Anti-Human Mrc2 Monoclonal Antibody For Flow Cytometry - Q9Ubg0

Item Number: CM0005736



Introduction

This -conjugated rabbit monoclonal antibody targets human MRC2 (CD280, ENDO180, UPARAP), a C-type mannose receptor involved in collagen internalization and matrix remodeling. Validated for flow cytometry applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human MRC2 mAb
Remarks/Alias	CD280; UPARAP; CLEC13E; ENDO180
Species	Human
GeneID (Human)	9902
GeneID	9902
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-513 of human MRC2 (NP_006030.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9UBG0
Protein Weight	167kDa
Shipping	Ice bag