

Rabbit Anti-Human Nectin-1/Cd111 Mab For Fc - Q15223

Item Number: CM0009579



Introduction

-conjugated rabbit monoclonal antibody targeting human Nectin-1/CD111, a cell adhesion molecule and herpesvirus entry receptor. Validated for flow cytometry, facilitating research in neurodevelopment and viral pathogenesis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Nectin-1/CD111 mAb
Remarks/Alias	ED4; PRR; HlgR; HV1S; HVEC; OFC7; PRR1; PVRR; CD111; PVRL1; PVRR1; SK-12; CLPED1; nectin-1
Species	Human
Gene ID (Human)	5818
Gene ID	5818
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-334 of human Nectin-1/CD111 (NP_002846.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q15223
Protein Weight	57kDa/51kDa/39kDa
Shipping	Ice bag