


Rabbit Anti-Human Nectin-1/Cd111 Monoclonal Antibody For Fc - Q15223

Item Number: CM0005282



Introduction

-conjugated rabbit monoclonal antibody against human Nectin-1/CD111 for flow cytometry. Validated for human samples, recognizes a 57/51/39 kDa protein. Suitable for cell adhesion, neurodevelopment, and virus entry research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Nectin-1/CD111 mAb
Remarks/Alias	ED4; PRR; HlgR; HV1S; HVEC; OFC7; PRR1; PVRR; CD111; PVRL1; PVRR1; SK-12; CLPED1; nectin-1
Species	Human
GeneID (Human)	5818
GeneID	5818
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-334 of human Nectin-1/CD111 (NP_002846.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q15223
Protein Weight	57kDa/51kDa/39kDa
Shipping	Ice bag