

Rabbit Anti-Human Nectin-4 Monoclonal Antibody For Fc - Q96Ny8

Item Number: CM0007741



Introduction

-conjugated rabbit monoclonal antibody targeting human Nectin-4 (PVRL4), validated for flow cytometry (FC). Suitable for research on cell adhesion, breast carcinoma, and measles virus receptor.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Nectin-4 mAb
Remarks/Alias	LNIR; PRR4; EDSS1; PVRL4; nectin-4
Species	Human
GeneID (Human)	81607
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-146 of human Nectin-4 (NP_112178.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q96NY8
Protein Weight	24kDa/55kDa
Shipping	Ice bag