

Rabbit Anti-Human Nt5E/Cd73 Monoclonal Antibody For Fc - P21589

Item Number: CM0012408



Introduction

conjugated rabbit monoclonal antibody against human NT5E/CD73, validated for flow cytometry (FC). Targets ecto-5'-nucleotidase, a GPI-anchored membrane enzyme hydrolyzing AMP. Suitable for immuno-oncology and nucleotide metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human NT5E/CD73 mAb
Remarks/Alias	NT; eN; NT5; NTE; eNT; CD73; E5NT; CALJA
Species	Human
Gene ID (Human)	4907
Gene ID	4907
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-547 of human NT5E/CD73 (NP_002517.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21589
Protein Weight	63kDa
Shipping	Ice bag