


Rabbit Anti-Human Opn1Lw/Opn1Mw Mab - P04000 / P04001

Item Number: CM0011623



Introduction

High-performance -conjugated rabbit monoclonal antibody targeting human OPN1LW/OPN1MW opsins for flow cytometry. Ideal for color vision and photoreceptor research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human OPN1LW/OPN1MW mAb
Remarks/Alias	CBP; RCP; ROP; CBBM; COD5; CBD; GCP; GOP; OPN1MW1
Species	Human
GeneID (Human)	5956/2652
GeneID	5956/2652
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-48 of human OPN1LW/OPN1MW (NP_064445.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P04000/P04001
Protein Weight	40kDa
Shipping	Ice bag