


Rabbit Anti-Human Opn1Lw/Opn1Mw Monoclonal Antibody For Flow Cytometry - P04000 / P04001

Item Number: CM0004115



Introduction

-conjugated rabbit monoclonal anti-human OPN1LW/OPN1MW antibody for flow cytometry. Targets red-sensitive (OPN1LW) and green-sensitive (OPN1MW) cone opsins, suitable for retinal research and color vision studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human OPN1LW/OPN1MW mAb
Remarks/Alias	CBP; RCP; ROP; CBBM; COD5; CBD; GCP; GOP; OPN1MW1
Species	Human
GeneID (human)	5956/2652
GeneID	5956/2652
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-48 of human OPN1LW/OPN1MW (NP_064445.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P04000/P04001
Protein Weight	40kDa
Shipping	Ice bag