


Rabbit Anti-Human P-Cadherin Mab - P22223

Item Number: CM0010429



Introduction

-conjugated rabbit monoclonal antibody against human P-cadherin (CDH3) for flow cytometry. Detects human P-cadherin in epithelial tissues and carcinomas, ideal for cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human P-cadherin mAb
Remarks/Alias	CDHP; HJMD; PCAD
Species	Human
Gene ID	1001
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 108-654 of human P-cadherin (NP_001784.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P22223
Protein Weight	87kDa/91kDa
Shipping	Ice bag