


Rabbit Anti-Human P-Cadherin Monoclonal Antibody For Flow Cytometry - P22223

Item Number: CM0009580



Introduction

-conjugated rabbit monoclonal antibody targeting human P-cadherin (CDH3) for flow cytometry. Recognizes cadherin-3 in epithelial tissues and carcinoma cells; useful for cell adhesion and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human P-cadherin mAb
Remarks/Alias	CDHP; HJMD; PCAD
Species	Human
GeneID	1001
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 108-654 of human P-cadherin (NP_001784.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P22223
Protein Weight	87kDa/91kDa
Shipping	Ice bag