

Rabbit Anti-Human P4Hb Monoclonal Antibody For Intracellular Flow Cytometry - P07237

Item Number: CM0009674



Introduction

-conjugated rabbit monoclonal antibody targeting human P4HB (protein disulfide-isomerase). Validated for intracellular flow cytometry, detecting a 57 kDa protein. Suitable for protein folding and ER stress studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human P4HB mAb
Remarks/Alias	DSI; GIT; PDI; PHDB; PDIA1; PO4DB; PO4HB; PROHB; CLCRP1; ERBA2L; P4Hbeta
Species	Human
Gene ID (Human)	5034
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 329-508 of human P4HB (NP_000909.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P07237
Protein Weight	57kDa
Shipping	Ice bag