

Rabbit Anti-Human Pcsk9 Monoclonal Antibody For Flow Cytometry - Q8Nbp7

Item Number: CM0012644



Introduction

-conjugated rabbit monoclonal antibody targeting human PCSK9 (Q8NBP7). Validated for intracellular flow cytometry applications. Recombinant immunogen corresponding to amino acids 31-692 of human PCSK9. Ideal for cholesterol metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human PCSK9 mAb
Remarks/Alias	FH3; PC9; FHCL3; NARC1; LDLCQ1; NARC-1; HCHOLA3
Species	Human
GeneID (Human)	255738
GeneID	255738
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 31-692 of human PCSK9 (NP_777596.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	Q8NBP7
Protein Weight	74kDa
Shipping	Ice bag