

Rabbit Anti-Human Perforin Mab For Fc (Intra) - P14222

Item Number: CM0009936

Introduction

-conjugated rabbit monoclonal antibody targeting human Perforin (PRF1) for intracellular flow cytometry. Validated for human samples; detects a 61 kDa protein. Suitable for immune cell cytotoxicity research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Perforin mAb
Remarks/Alias	P1; PFP; HPLH2
Species	Human
Gene ID (Human)	5551
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-555 of human Perforin (NP_005032.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P14222
Protein Weight	61 kDa
Shipping	Ice bag