


Rabbit Anti-Human Phospho-PI3K P85A (Y467)/P85B (Y464)/P55γ (Y199) Monoclonal Antibody For Flow Cytometry (Intracellular) - P27986 / O00459 / Q92569

Item Number: CM0001539



Introduction

-conjugated rabbit monoclonal antibody specific for phosphorylated PI3K regulatory subunits p85α (Y467), p85β (Y464), and p55γ (Y199). Validated for intracellular flow cytometry to analyze PI3K/AKT signaling in human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Phospho-PI3K P85α-Y467+P85β-Y464+P55γ-Y199 mAb
Remarks/Alias	PIK3R1/PIK3R2/PIK3R3
Species	Human
Gene ID (Human)	5295/5296/8503
Gene ID	5295/5296/8503
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P27986/O00459/Q92569
Protein Weight	42-54kDa/81-84kDa
Shipping	Ice bag