

Rabbit Anti-Human Phospho-Pi3Kp85A-Y467+P85B-Y464+P55γ-Y199 Mab - P27986 / O00459 / Q92569

Item Number: CM0001913



Introduction

-conjugated rabbit monoclonal antibody against human phospho-PI3K p85α (Y467), p85β (Y464) and p55γ (Y199). Validated for intracellular flow cytometry (FC). Ideal for studying PI3K/AKT signaling in human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Phospho-PI3Kp85α-Y467+P85β-Y464+P55γ-Y199 mAb
Remarks/Alias	PIK3R1/PIK3R2/PIK3R3
Species	Human
GeneID (Human)	5295/5296/8503
GeneID	5295/5296/8503
Immunogen	Synthetic peptide / A synthetic phosphorylated peptide around Y467 of human PI3Kp85α (NP_852664.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P27986/O00459/Q92569
Protein Weight	42-54kDa/81-84kDa
Shipping	Ice bag