

Rabbit Anti-Human Plzf Monoclonal Antibody For Intracellular Flow Cytometry - Q05516

Item Number: CM0012302



Introduction

-conjugated rabbit monoclonal antibody targeting human PLZF (ZBTB16), a transcriptional repressor involved in myeloid maturation. Validated for intracellular flow cytometry.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human PLZF mAb
Remarks/Alias	PLZF; ZNF145
Species	Human
Gene ID (Human)	7704
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 6-137 of human PLZF (NP_005997.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	Q05516
Protein Weight	74kDa
Shipping	Ice bag