

## Rabbit Anti-Human Plzf Monoclonal Antibody For Intracellular Flow Cytometry - Q05516

Item Number: CM0009307



### Introduction

-conjugated rabbit monoclonal antibody targeting human PLZF (ZBTB16). Validated for intracellular flow cytometry. Reactivity: Human. Explore transcriptional repression and myeloid differentiation mechanisms.

[Learn More](#)

| Parameter        | Specification                                                         |
|------------------|-----------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human PLZF mAb                                            |
| Remarks/Alias    | PLZF; ZNF145                                                          |
| Species          | Human                                                                 |
| Gene ID          | 7704                                                                  |
| Immunogen        | Recombinant fusion protein containing amino acids 6-137 of human PLZF |
| Source           | Rabbit                                                                |
| Category         | Monoclonal Antibodies                                                 |
| Application      | FC (intra)                                                            |
| Cross Reactivity | Human                                                                 |
| SWISS            | Q05516                                                                |
| Protein Weight   | 74kDa                                                                 |
| Shipping         | Ice bag                                                               |