

Rabbit Anti-Human Podoplanin Monoclonal Antibody For Flow Cytometry - Q86YI7

Item Number: CM0012330



Introduction

-conjugated rabbit monoclonal antibody targeting human Podoplanin (PDPN/Aggrus/Gp36). Validated for flow cytometry, human cross-reactivity. Suitable for cancer biology, lymphatics, and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Podoplanin mAb
Remarks/Alias	PDPN; AGGRUS; GP36; GP40; Gp38; HT1A-1; OTS8; PA2.26; T1A; T1A-2; T1A2; TI1A; podoplanin
Species	Human
GeneID (Human)	10630
GeneID	10630
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-123 of Human Podoplanin (NM_006474.5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q86YL7
Protein Weight	12kDa/16kDa/18kDa/24kDa
Shipping	Ice bag