


Rabbit Anti-Human Pth1R Monoclonal Antibody () For Flow Cytometry - Q03431

Item Number: CM0010046



Introduction

Rabbit monoclonal antibody against human PTH1R, conjugated to , validated for flow cytometry. Parathyroid hormone 1 receptor (PTH1R) is a G protein-coupled receptor involved in calcium homeostasis. Ideal for research and diagnostic applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human PTH1R mAb
Remarks/Alias	PFE; EKNS; PTHR; PTHR1
Species	Human
Gene ID	5745
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-188 of human PTH1R (NP_000307.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q03431
Protein Weight	66kDa
Shipping	Ice bag